

# Everlasting ABI 310

**Extend the life span of your single capillary sequencer**  
**ABI Prism® 310**



**We offer a full service contract, Win10 upgrade included!**

## OUR OFFER INCLUDES

**(THROUGHOUT THE DURATION OF THE CONTRACT)**

- » Inspection / Preventive examination of your sequencer yearly, including operational and process qualification (OQ/PQ) and instrument performance verification (IPV)
- » Defect diagnosis and all repairs
- » New data analysis computer: Windows 10 license, installation of Data Collection and analysis software
- » Labour, parts and travel costs of service engineer
- » Unlimited application support

**NO NEED TO BUY ADDITIONAL SOFTWARE LICENSES, KEEP ON USING YOUR CURRENT ONES!**

### Please notice:

The offer is valid for ABI Prism® 310 sequencer only, co-installed with a data collection computer with Windows NT Workstation, Windows 2000 or Windows XP and includes delivery of the new data analysis PC with Windows 10. Installation of data collection and analysis software (Data Collection v3.X, GeneScan from v3.X, Sequencing Analysis from v3.X, GeneMapper from v3.X/ID) is also included. This new PC is put into operation during first regular inspection of the sequencer and is covered by warranty throughout the duration of the contract. Apart from Windows 10 this offer does not cover any additional software licenses.