

# Upgrading ABI 3130/3130xl

**Extend  
the life span  
of your  
sequencers  
ABI 3130  
and 3130xl!**



**Windows 10 upgrade,  
standalone or as a part  
of service contract.**

## OUR OFFER INCLUDES:

- » New data analysis computer: Windows 10 license
- » Installation of Data Collection and analysis software
- » Preventive maintenance of your instrument including the operational and process qualification (OQ/PQ) and instrument performance verification (IPV)
- » Labour and travel costs of service engineer
- » The upgrade can be performed as a part of long-term service contract

**NO NEED TO BUY ADDITIONAL SOFTWARE LICENSES,  
KEEP ON USING YOUR CURRENT ONES!**

The offer is valid for ABI 3130/3130xl sequencers only, co-installed with a data collection computer with Windows XP or older OS and includes delivery of the new data analysis PC with Windows 10 and installation of data collection and analysis software (DataCollection and any version of Sequencing Analysis, GeneMapper, SeqScape / Variant Reporter). This new PC is put into operation during first regular inspection and is covered by 2-years warranty. Apart from Windows 10 this offer does not cover any additional software licenses.