

AFTER OUR MAINTENANCE YOUR SEQUENCERS BLOSSOM

Instrument validation of Applied Biosystems® sequencers, all models

Performed according to "Best practice guidelines" of CMGS (Clinical Molecular Genetics Society)

Operational and Process Qualification (OQ/PQ)

- » The OQ/PQ Protocol verifies that the Applied Biosystems genetic analyzer, at the time of testing, is operating in accordance with its specifications.

DNA Sequencer Process Validation and Verification (IPV)

- » The Instrument Performance Verification procedure demonstrates expected performance of the instrument has been met, effectively and reproducibly.
- » The verification is performed for DNA sequencing and fragment analysis, respectively, using validation standards of the manufacturer.

**At the regular inspection
the OQ/PQ/IPV steps are
performed**

FOR FREE!

Applied Biosystems is a registered trademark of Life Technologies.

Order your instrument inspection
including validation at service@seqme.eu.

SEQme s. r. o., Dlouhá 176, Dobruška, Czech republic

service@seqme.eu

SEQme s.r.o. is an independent consulting and training company and a service provider in the field of DNA sequencing and Real-Time PCR.

More information about our services, trainings in Real-Time PCR, DNA sequencing and other techniques can be obtained at www.seqme.eu.

© 2015 SEQme s. r. o. | All rights reserved.

 **SEQme**
SEQ & QPCR
COMPANY